

URINE THC

The EMIT II Plus Cannabinoid Assay detects the major metabolite of Δ^9 -THC, 11-nor- Δ^9 -THC-9-carboxylic acid. It also detects other Δ^9 -THC metabolites.

Non Interfering Substances

Each of the following compounds when added to urine at +/- 25% concentration of the cutoff do not yield a false response relative to the 50 ng/mL cutoff:

Compound	Concentration
Acetone	1.0 g/dL
Ascorbic acid	1.5 g/dL
Bilirubin	0.25 mg/dL
Creatinine	0.5 g/dL
Ethanol	1.0 g/dL
Gamma Globulin	0.5 g/dL
Glucose	2.0 g/dL
Hemoglobin	115 mg/dL
Human serum albumin	0.5 g/dL
Oxalic Acid	0.05 g/dL
Riboflavin	7.5 mg/dL
Sodium Chloride	1.5 g/dL
Urea	6.0 g/dL

Specificity

The table below lists the concentration of compounds that produce results approximately equivalent to the 50 ng/mL cutoff. Each concentration represents the reactivity level for the stated compound when the compound is added to a negative urine specimen. These concentrations are within the range of the levels found in urine following use of the drug or, in the case of metabolites, the parent compound. If a specimen contains more than one compound detected by the assay, lower concentrations than those listed may combine to produce a rate approximately equivalent to or greater than the cutoff calibrator of 50 ng/mL.

Compound	20 ng/mL cutoff	50 ng/mL cutoff	100 ng/mL cutoff
8- β -11-Dihydroxy- Δ^9 -THC	24	58	109
8- β -11-Hydroxy- Δ^9 -THC	26	68	146
11- Hydroxy- Δ^8 -THC	43	67	129
11- Hydroxy- Δ^9 -THC	42	77	124
9-Carboxy-11-nor- Δ^9 -THC-glucuronide	79	95	328
Δ^8 -THC	79	220	660
Δ^9 -THC	78	220	620

The table below lists the concentrations of compounds that show a negative response to the Emit II Plus THC Assay at all cutoff levels. Positive results for specimens containing other compounds structurally unrelated to cannabinoids have not been observed.

Compound	Concentration
Acetaminophen	1000 µg/mL
α-Acetyl-N,N-dinormethadol (dinor LAAM)	25 µg/mL
L-α-Acetylmethadol (LAAM)	25 µg/mL
N-Acetylprocainamide (NAPA)	400 µg/mL
Acetylsalicylic Acid	1000 µg/mL
Amitriptyline	1000 µg/mL
D-Amphetamine	1000 µg/mL
Benzoylcegonine	1000 µg/mL
Buprenorphine	100 µg/mL
Caffeine	1000 µg/mL
Cimetidine	1000 µg/mL
Clomipramine	2.5 µg/mL
Clonidine	1000 µg/mL
Codeine	500 µg/mL
Cotinine	100 µg/mL
Cyclobenzaprine	1000 µg/mL
Desipramine	800 µg/mL
Diphenhydramine	1000 µg/mL
Doxepin	1000 µg/mL
2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine (EDDP)	1000 µg/mL
Fluoxetine	1000 µg/mL
Glutethimide	500 µg/mL
Ibuprofen	1000
Ketamine	100 µg/mL
Ketorolac Tromethamine	1000 µg/mL
Lormetazepam	1 µg/mL
LSD	10 ng/mL
Meperidine	1000 µg/mL
D-Methamphetamine	35 µg/mL
Methaqualone	1500 µg/mL
Morphine	1000 µg/mL
Naproxen	1000 µg/mL
Nortriptyline	1000 µg/mL
Oxazepam	300 µg/mL
Phencyclidine	1000 µg/mL
Phenytoin	1000 µg/mL
Promethazine	1000 µg/mL
Propoxyphene	1000 µg/mL
Ranitidine	1000 µg/mL
Scopolamine	500 µg/mL
Secobarbital	1000 µg/mL
Thioridazine	100 µg/mL
Tramadol	1000 µg/mL

Tyramine	100 µg/mL
Zidovudine (AZT)	2 mg/mL
Zolpidem	100 µg/mL

Sensitivity

The sensitivity level (minimum detection limit) of the Emit II Plus Cannabinoid Assay using the 50 ng/mL cutoff is 35 ng/mL. This level represents the lowest level of Δ^9 -THC, 11-nor- Δ^9 -THC-9-carboxylic acid that can be distinguished from 0 ng/mL with a confidence level of 95%.

CAMC laboratories use the 50 ng/mL cutoff for the cannabinoid assay.

(Syva/EMIT are trademarks of **Siemens Healthcare** Diagnostics Inc.)

